



Apply the Update

- 1) Once you have saved the update file on your computer, go to the location of the file. Downloaded files are commonly saved in your **My Downloads** folder, by default. Once you have found the downloaded **ABAS3.exe** file, right-click on it and choose **Copy**. Open the **C:\WPS\ABAS31** folder. Right-click anywhere inside this folder and choose **Paste**.



- 2) Choose **Copy and Replace**, when prompted by Windows.



- 3) Proceed to plug in your **ABAS-3 USB** Key and wait for a solid green light on the tip of it. Continue to open up your **ABAS-3** program.